

stethoscopes

Stethoscopes	429 - 436
Baby/Neonatal Stethoscopes	430 - 431
Electronic Stethoscopes	432
Classic Stethoscopes	433 - 434
Cardiology Stethoscopes	435

ALL PRICES
EXCLUDE VAT

info@henryscheinmedical.co.uk

phone 0800 0304169

429

Stethoscopes

Tristar

This high quality stethoscope has three different sized double chest-pieces for adult, children, babies, and new born and is suitable for perfect auscultation in all frequency ranges.



LATEX FREE

Riester
The familiar way



Neon



Baby



Duplex

- Double chest pieces for Duplex, Duplex Baby, and Duplex Neonatal, made of aluminium with rapid exchange connector
- Special membranes diameter 48mm, 36mm, 28mm for first class acoustics
- Smooth adaption of the chest piece with the help of non-chill rims
- Bells diameter 36mm, 28mm, 24mm
- Matt chromium plated binaural with integrated reinforced double spring
- Optimal sealing of the external auditory canal with soft, replaceable eartips with metal screw thread
- Strong Y tube made of PVC for optimal acoustic to shield against undesired stray sound
- Tristar is supplied with two pairs of replacement eartips a replacement membrane and a name plate
- Grey only

STE250

£62.70

Duplex de Luxe



LATEX FREE

Riester
The familiar way

- Double chest piece made of either stainless steel or light weight aluminium
- Special membranes diameter 48mm and 36mm
- Smooth adaption of the chest piece with the help of non-chill rims
- Strong Y tube made of PVC for optimal acoustic to shield against undesired stray sound
- Overall length: 77cm
- Matt chromium plated binaural with integrated reinforced double spring
- Extra soft, replaceable eartips with metal screw thread for comfortable and effective sealing of the external auditory canal
- Supplied with two pairs of replacement eartips a replacement membrane and a name plate
- Choice of three colours Slate Grey, Saffron, and Green

STE251

Specify Colour

£45.60

Spirit Cardiology

5 years
Guarantee



- High quality stainless steel combination chestpiece for excellent acoustics
- Dual head
- 635mm Tubing
- 2 spare pairs of soft seal eartips and 1 hard
- Black only

STE037

£96.00

Spirit Classic

3 years
Guarantee



- High quality stainless steel chestpiece
- Dual head
- 710mm Tubing
- 2 spare pairs of soft seal eartips and 1 hard
- Available in Blue, Black, Grey or Burgundy
- One piece binaural with adjustable headframe and integral spring

STE035

Adult - specify colour

£32.50

STE036

Paediatric - specify colour

£38.00

Paediatric Stethoscope



• Colour Black



STE058 | Paediatric Stethoscope | **£9.50**

Maxima Stethoscope Dualhead Adult



• Choice of Burgundy, Black or Navy Blue colours



STE054 | Maxima Stethoscope Dualhead Adult, Specify colour | **£29.95**

Infant Stethoscope



• Colour Grey



STE059 | Infant Stethoscope | **£9.50**

Maxima Stethoscope Rappaport




• Choice of Black, Navy Blue, Neon Orange, or Royal Blue colours



STE055 | Maxima Stethoscope Rappaport, Specify colour | **£12.50**

Dualhead Stethoscope



• Choice of Black, grey, Royal Blue and Red colours



STE057 | Dualhead Stethoscope, Specify colour | **£9.50**

Nurse Stethoscope



• Choice of Black, grey, Royal Blue and Red colours



STE056 | Nurse Stethoscope, Specify colour | **£7.50**

Stethoscopes

Littmann® Electronic 3100 Stethoscope



- 85% Ambient Noise Reduction
- 24 x sound amplification
- Looks and feels like a traditional stethoscope
- LCD user interface with easy to operate controls, battery level indicator and heart rate monitor
- Operates on a single AA battery
- Auto-on feature; unit stays in battery-saving 'sleep mode'
- 15mm diameter sound sensor suitable for infants
- The solid metal chestpiece and rugged design offers years of trouble-free operation
- Clinical evidence shows that it's easier to detect difficult-to-hear heart sounds, like S3 gallops, aortic regurgitation murmurs, and abnormal lung sounds

STE4220	Black	£259.00
STE42202	Burgundy	£259.00
STE42203	Navy Blue	£259.00
STE42204	Hunter Green	£259.00
STE42205	Plum	£259.00
DIS352	AA Batteries, pack of 4	£2.97



Littmann® Electronic 3200 Stethoscope



- Same as 3100 Model with additional features
- FREE Bluetooth® dongle included
- FREE Zargis® StethAssist™ Heart and Lung Sound Visualization Software
- Record & playback function in bell, diaphragm or extended range
- Display and print sounds as a phonocardiogram
- Send as an e-mail attachment or connect to optional ZargisTelemed.com to share recordings online with colleagues across the globe
- Facilitate storage of sounds in patient records or transfer to colleagues for further analysis
- 2 years warranty
- 69cm - 27"

STE42301	Black	£320.00
STE42302	Burgundy	£320.00
STE42303	Navy Blue	£320.00
DIS352	AA Batteries, pack of 4	£2.97



Littmann® Classic II SE

- Traditional two-sided chestpiece with bell/diaphragm
- Tuneable diaphragm for high and low frequencies
- High performance single lumen tubing
- Stainless steel chestpiece
- Non-chill rim
- 3 years warranty
- 71cm - 28"

3M

3 years
Guarantee

LATEX FREE



STE110	71cm, Black tubing	£55.00
STE100	71cm, Orange tubing 1	£57.50
STE104	71cm, Grey tubing	£57.50
STE103	71cm, Navy Blue tubing	£57.50
STE108	71cm, Ceil Blue tubing 2	£57.50
STE109	71cm, Seafoam Green tubing 3	£57.50
STE119	71cm, Hunter Green tubing	£57.50
STE102	71cm, Burgundy tubing	£57.50
STE120	71cm, Caribbean Blue tubing	£57.50
STE139	71cm, Baby Pink tubing	£57.50
STE111	71cm, Purple tubing	£57.50
STE112	71cm, Rainbow tubing	£57.50
STE118	71cm, Brass tubing	£57.50
STE121	71cm, Copper tubing	£57.50
STE806	Name Tag, Black	£5.65



Paediatric Version

• 25mm Bell

STE106	71cm, Black tubing	£66.00
STE105	71cm, Red tubing	£66.00
STE113	71cm, Caribbean Blue tubing	£66.00
STE137	71cm, Raspberry tubing 4	£66.00
STE114	71cm, Rainbow tubing	£66.00
STE806	Name Tag, Black	£5.65


Infant Version

• 19mm Bell

STE107	71cm, Black tubing	£66.00
STE116	71cm, Red tubing	£66.00
STE138	71cm, Raspberry tubing	£66.00
STE117	71cm, Lilac tubing 5	£66.00
STE115	71cm, Rainbow tubing	£66.00
STE806	Name Tag, Black	£5.65



Name Tag
Black identification tag includes one written-on and one aluminium name plate for engraving



Littmann® Select



LATEX FREE

3M

2 years
Guarantee

- Single sided chestpiece
- Tuneable diaphragm for high and low frequencies
- Ideal stethoscope for nurses and technicians
- 2 years warranty
- 71cm - 28"




STE162	71cm, Burgundy tubing	£33.00
STE163	71cm, Navy Blue tubing	£33.00
STE165	71cm, Caribbean, Blue tubing	£33.00
STE166	71cm, Lilac tubing	£33.00
STE806	Name Tag, Black	£5.65

Littmann® Teaching Stethoscopes



3M



- The pupil can share the sounds with a teacher at the same time with the convenience of a double headset
- Two-sided chestpiece
- Tuneable diaphragm on one side augments the traditional bell function of the opposite side

- Stainless steel chestpieces
- High performance single-lumen tubing
- Non-chill rim
- Black only

STE170		£89.00
--------	--	--------

Stethoscopes

Littmann® Master Classic II

LATEX FREE

3M



3 years
Guarantee

Master Classic II all Black edition (STE129) Black plated chestpiece and eartips

- Tuneable diaphragm for high and low frequencies
- Single sided chestpiece
- High performance single lumen tubing
- Non-chill rim
- Comfortable, durable and portable
- 3 years warranty
- 69cm - 27"

STE123	69cm, Black tubing	£75.00
STE124	69cm, Burgundy tubing	£75.00
STE125	69cm, Navy Blue tubing	£75.00
STE128	69cm, Caribbean Blue tubing	£75.00
STE129	69cm, All Black edition	£75.00
STE806	Name Tag, Black	£5.65

NEW Littmann® Lightweight Nurse Stethoscope

3M

The choice for nurses

This general-purpose, entry-level stethoscope is an excellent, affordable choice for blood-pressure and patient assessment. Littmann quality is apparent from the soft-sealing eartips to the combination bell/diaphragm chestpiece.

- Provides clear, reliable acoustic performance
- Ovoid chestpiece with offset bell for better orientation around blood-pressure cuffs and body contours
- Tunable diaphragm for easy low- and high-frequency monitoring
- Adjustable and anatomically correct headset
- Comfortable to use
- 3M Littmann® Snap Tight Soft-Sealing Eartips for an excellent acoustic seal
- Comfortable to use
- Lightweight and durable design
- Choice of colours - Ceil Blue, Burgundy & Black



STE167	L/Weight II SE Stethoscope Nursing Black	£24.95
STE168	L/Weight II SE Stethoscope Nursing C/Blue	£24.95
STE169	L/Weight II SE Stethoscope Nursing Burg	£24.95

Littmann® Master Cardiology Stethoscope



3M



Master Cardiology all Black edition (STE150)
Black plated chestpiece and eartubes

7 years
Guarantee

LATEX FREE

- Handcrafted, solid stainless steel chestpiece
- Outstanding acoustic quality
- Single-sided tuneable diaphragm for high and low frequencies
- Double lumen tubing
- Special procedures adaptor for difficult site and paediatric auscultation
- Non-chill rim
- 7 years warranty
- 69cm - 27"

STE145	69cm, Black tubing	£160.00
STE146	69cm, Burgundy tubing	£160.00
STE151	69cm, Plum tubing	£160.00
STE150	69cm, All Black edition	£160.00
STE806	Name Tag, Black	£5.65

Littmann® Cardiology S.T.C. Stethoscope



3M

6 years
Guarantee

LATEX FREE



- Unique patented chestpiece shape
- Outstanding acoustic quality
- Easy to grasp
- Ovoid chestpiece, allows easy movement between auscultation sites
- Single-sided tuneable diaphragm for high and low frequencies
- Double lumen tubing
- 6 years warranty
- 69cm - 27"

STE154	69cm, Black tubing	£155.00
STE155	69cm, Navy Blue tubing	£155.00
STE156	69cm, Burgundy tubing	£155.00
STE157	69cm, Plum tubing	£155.00
STE806	Name Tag, Black	£5.65

Littmann® Cardiology III Stethoscope



3M

5 years
Guarantee

LATEX FREE




- Combination chestpiece
- Large side for adult patients
- Small side for paediatric patients
- Tuneable diaphragm on each side of the chestpiece
- Possible conversion of small side into a traditional bell by removing the diaphragm and replacing with the non-chill rim
- Double lumen tubing
- 5 years warranty
- 69cm - 27"

STE132	69cm, Black tubing	£138.00
STE133	69cm, Burgundy tubing	£138.00
STE134	69cm, Navy Blue tubing	£138.00
STE136	69cm, Plum tubing	£138.00
STE131	69cm, Brass Edition	£138.00
STE140	69cm, Copper Edition	£138.00
STE806	Name Tag, Black	£5.65

Equipment Financing

HENRY SCHEIN®
FINANCIAL SERVICES

Finance Solutions for Healthcare Practitioners



Build a Secure Future with Practice Financing Options from Henry Schein Financial Services

Finance Solutions for Healthcare Practitioners

Our low-cost, full service equipment and practice financing capabilities give you the options you need to move your business forward. Our nationwide specialists offer their expert advice and support with your best interests in mind.

Low Rates

Pre-Approval up to
£75,000 with no financials*

24 hour approvals

Flexible programmes

Terms from 1-7 years

First Class Service

Free up existing credit lines
and working capital

Lease, Hire Purchase
and Business Loans

de lage landen 
partners in finance
Preferred supplier of the NHS



* Finance offered is subject to status and credit approval

Get your personalised quote today...

01634 877540

HENRY SCHEIN®
FINANCIAL SERVICES in partnership with

ChoiceHealth
LEASING

www.henryschein.co.uk